

# Eight characteristics Concern about your need



Display blood pressure



One click measurement



HD screen



2\*99 groups Memory



Self-measure with cuff



Voice broadcast



Remind arrhythmia



Divide into A/B groups

### **Product Details**

Model	HF-A21	
Operation	Temperature:5~40°C/41-104°F	
	Humidity:15%~85%RH	
Storage	Temperature:-20°C~55°C / -4 ~ 131°F	
	Humidity:10%-95%RH	
N.W	About 120 g	
Product	67(L)×66(W)×28(H)mm	
dimensions		
Measuring Range	Blood Pressure:0 to 280 mmHg( 0 to 37.3 kPa)	
	Pulse Rate:40 to 195 times/minute	
Service life	5 Years (Without Battery)	
Display Mode	Digital liquid crystal display	
Accuracy	Blood Pressure: ± 3mmHg /±0.4kpa	
	Pulse: ±5%	
Power Supply	2*AAA Alkaline battery	
Enclosure	Monitor*1,certificate*1,Battery*2,Stock case*1,user manual*1	



## **One Key Operation**






## 2\*99 groups memory

2 \* 99 groups , automatic storage of measurement resulets , more convenience than manual storage







### **IHB HEART RATE TEST**

In the measurement of blood pressure at the same time one the heart rate is also tested, a machine is with multi-function





## SMART VOICE BROADCAST (OPTION)-SUPPORT YOUR LOCAL LANGUAGE

Support your local Language, for the Elderly with known date difficulties. voice broadcast can convey blood pressure information more conveniently






### STOCK CASE

Put in case and easy to carry to outside for test




# BLOOD PRESSURE CLASSIFICATION (WHO) INDICATOR

Blood pressure classification (WHO) indicator






# YOU CAN DO IT AS NEEDED TURN SOUND ON OR OFF





## AUTOMATIC POWER OFF AND LOW BATTERY REMINDE



#### OPERATING INSTRUCTIONS

Improve the accuracy of star measurement-wrist-type electronic blood pressure monitor correctly measures the pose code

The correct measurement posture is also very important to the accuracy of the measurement results

- Straighten your body and maintain a correct sitting posture.
- \*Be careful not to bend your body forward
- Take a deep breath 5-6 times before the measurement, restore natural breathing, relax the shoulders and arms, and relax the whole body
- When measuring, keep the wrist strap consistent with the height of the heart (the position of the nipple)

The position sensing function will tell you the ideal measuring height.



#### Correct wrist posture

Don't force your fingers, straighten them naturally




Improving the measurement of Philip-wrist type electronic sphygmomanometer wrong measurement posture

Please aoid the following postures or gestures when measuring



- Even if the position of the wristband is correct, the following actions may still cause the measurement value to be higher, especially when you are holding the head tightly
- Do not support the wristband with your other hand, otherwise it will cause false detection



#### **CERTIFICATE**











#### **COMPANY PROFILE**

Welcome to our company





Shenzhen Hopfat Medical Instrument Co. Ltd was established in 2010, located in Baoan District ,Shenzhen, China. Hopfat is a professional global medical device OEM/ODM manufacturing supplier in researching, manufacturing and exporting which mainly engaged in digital electronic Blood pressure monitor, Digital electronic thermometers, Non-contact forehead thermometer and other medical devices. The company has passed CE, RHOS, FCC and other certifications. At the same time; it has more than 10 patents and a number of software copyrights; the factory strictly in accordance with ISO13485.

quality management system for production management, the current products have been exported to more than 50 countries and region. Such as United States, Netherlands, Kenya, Vietnam.

Hopfat has a professional design team who make design as customer's requirement, core software algorithms and clinical improvement capabilities; it also has completely automatic equipment and multiple assembly lines, the annual maximum production capacity of electronic blood pressure monitor is 2.4 million units, the infrared thermometers is 72 thousand units, and the Non-contact forehead thermometer is 1.2 million units;

Fast shipment. Will ship the goods in 3-5 days under 500 pcs. Will ship the goods in about 2 weeks

About After-service, Support 1-2 years quality guarantee in difference goods. We will be patient to answer or solve the questions that you have.

Trust and Choose us! Take care your health and contact us!

#### **OUR SERVICES**



#### OEM/ODM

1 OEM your native language. Such as Japanese, Korean, German, Spanish, Russian and so on.

2 Embossing your logo. Send your logo file to us, then we will make desing for your reference and then start to producdtion after confirm.





#### Order goods with our brand

Support to send some promotion or attractive pictures to you for publicity

#### **OUR TEAM**

Product R & D department









### **OUR PRODUCTION LINE**

Product quality inspection and assembly department







## **QUALITY CHECK SCENARIO**







the accuracy



## **PACKAGE DETAILS**







Package size:420\*250\*340mm Product Model:HF-A21 Quantity:24pcs/carton

GW:8KG




#### TRANSPORTATION DETAILS











Ship by our agent (DAP/CIF term):

- 1.Ship out the goods,
- Send the Tracking number and Tracking link to you,
- 3.Remind you when you are about to receive the goods.

Ship by your agent in China (FCA term):

- 1.Ship out the goods,
- Send the Tracking number and Tracking link to you
- 3.Remind you that your Chinese forwarder have already received the goods.

## **FAQ**

- 1. Are you Manufacturer or Trade company?
- ---We are a Manufacturer. You can check our factory video on Alibaba or you can video call with us to visit our factory.

#### 2.Why choose us?

We are factory, can provide reasonable price and good after-service. Provide OEM&ODM service, Such as Customized your native language, Embossing your logo on items or change the outlook as your requirement.

- 2.What's certificate you have?
- --- CE, FCC, RHOS and ISO13485.
- 3.What's your MOQ? Could I order some samples at first time?
- ---The MOQ is 24pcs for Blood pressure monitor. The MOQ is 240pcs for digital thermometer.

As our support, you can order some samples first to check the quality and market if you want.

- 4. What's your payment terms?
- --- Accept TT(Bank transfer), Trade Assurance or others.
- 5.What's your quality guarantee period?
- --- 1 year from shipment.
- 6.What's the sample lead time&production lead time?
- ---Will ship out in 2 days for samples. The production time is about 5-15days for bulk order.
- 7.Do you support OEM/ODM order?
- ---Sure. Just contact or call us. Then we will make the goods as your requirement.

#### **CONTACT US**





+86 18124768286



Sales6@hopfat-w.com



www.lianuoyhealth.com



