

## HF-E601

**Forehead Thermometer** 



## Six characteristics Concern about your need



32 groups of memory values



High-definition large-screen display



one-key temperature measurement



Backlight display



Automatic shut-down



High temperature warning

#### **Product Details**

Product name	Digital infrared thermometer
Model	HF-E601
Colour	Green
Display	LCD Digital Display
Memory	32 Sets
Accuracy	±0.2°C
Measurement range	32.0~42.9°C
Measuring distance	1~3CM
Battery	2*AAA battery
Measurement time	≦2 s
Low voltage prompt	2.6 (±0.2) V
Ambient temperature	16°C∼35°C
Measuring method	Non-Contact Measurement



# LARGE BUTTON AND ONE BUTTON MEASURING, EASY TO USE!

• • •





### **32 SETS MEMORY**





### **ACCURE MEASUREMENT**

2 Seconds measurement







# 1-3CM NON-CONTA & SAFETY MEASUREMENT

NON CONTACT TEMPERATURE MEASUREMENT





#### **FEVER ALARM**

Automatic fever alarm when Temperature is high





#### SMART VOICE BROADCAST(OPTION)-SUPPORT YOUR LOCAL LANGUAGE

...





# Large LCD screen with backlight (clear and easy to read)

You can also easily read even at night





#### **SWITCH BETWEEN°C &°F UNIT**





# AUTOMATIC POWER OFF AND LOW BATTERY REMINDER

2\*AAA BATTERIES ARE USED FOR POWER SUPPLY



2\*AAA BATTERIES






## AUTOMATIC POWER OFF AND LOW BATTERY REMINDER



X

Power off automatically in 30 seconds. No action required



The power icon will be displayed in the Displayed when the device is low

#### **CERTIFICATE**



Still more

Defentita lim

量()量







#### **COMPANY PROFILE**

Welcome to our company





Shenzhen Hopfat Medical Instrument Co. Ltd was established in 2010, located in Baoan District ,Shenzhen, China. Hopfat is a professional global medical device OEM/ODM manufacturing supplier in researching, manufacturing and exporting which mainly engaged in digital electronic Blood pressure monitor , Digital electronic thermometers, Non-contact forehead thermometer and other medical devices. The company has passed CE, RHOS,FCC and other certifications.

At the same time; it has more than 10 patents and a number of software copyrights; the factory strictly in accordance with ISO13485 quality management system for production management, the current products have been exported to more than 50 countries and region. Such as United States, Netherlands, Kenya, Vietnam.

Hopfat has a professional design team who make design as customer's requirement, core software algorithms and clinical improvement capabilities; it also has completely automatic equipment and multiple assembly lines, the annual maximum production capacity of electronic blood pressure monitor is 2.4 million units, the infrared thermometers is 72 thousand units, and the Non-contact forehead thermometer is 1.2 million units; Fast shipment. Will ship the goods in 3-5 days under 500 pcs. Will ship the goods in about 2 weeks

About After-service, Support 1-2 years quality guarantee in difference goods. We will be patient to answer or solve the questions that you have.

Trust and Choose us! Take care your health and contact us!


CLICK TO LEARN >>>

#### **OUR SERVICES**



#### OEM/ODM

- 1 OEM your native language. Such as Japanese, Korean, German, Spanish, Russian and so on.
- 2 Embossing your logo. Send your logo file to us, then we will make desing for your reference and then start to producdtion after confirm.



#### Order goods with our brand

Support to send some promotion or attractive pictures to you for publicity

### **OUR TEAM**

Product R & D department









### **OUR PRODUCTION LINE**

Product quality inspection and assembly department







### **QUALITY CHECK SCENARIO**







#### **Product inspection process**



We will check the accuracy with Mercury thermometer before shipment.



High precision constant temperature groove.



### **PACKAGE DETAILS**







Package size: 600\*250\*400mm Product Model: HF-E601

Quantity:48pcs/carton

GW:9KG




#### TRANSPORTATION DETAILS











Ship by our agent (DAP/CIF term):

- 1. Ship out the goods,
- 2. Send the Tracking number and Tracking link to you,
- 3.Remind you when you are about to receive the goods.

Ship by your agent in China (FCA term):

- 1. Ship out the goods,
- Send the Tracking number and Tracking link to you
- Remind you that your Chinese forwarder have already received the goods.

### **FAQ**

- 1. Are you Manufacturer or Trade company?
- ---We are a Manufacturer. You can check our factory video on Alibaba or you can video call with us to visit our factory.

#### 2.Why choose us?

We are factory, can provide reasonable price and good after-service. Provide OEM&ODM service, Such as Customized your native language, Embossing your logo on items or change the outlook as your requirement.

- 2.What's certificate you have?
- --- CE, FCC, RHOS and ISO13485.
- 3.What's your MOQ? Could I order some samples at first time?
- ---The MOQ is 24pcs for Blood pressure monitor. The MOQ is 240pcs for digital thermometer.

As our support, you can order some samples first to check the quality and market if you want.

- 4. What's your payment terms?
- --- Accept TT(Bank transfer), Trade Assurance or others.
- 5.What's your quality guarantee period?
- --- 1 year from shipment.
- 6.What's the sample lead time&production lead time?
- ---Will ship out in 2 days for samples. The production time is about 5-15days for bulk order.
- 7.Do you support OEM/ODM order?
- ---Sure. Just contact or call us. Then we will make the goods as your requirement.

#### **CONTACT US**





+86 18124768286



Sales6@hopfat-w.com



www.lianuoyhealth.com



